CELESTRON[®]
LABS
DIGITAL
ADVANCED

LCD Desktop Microscope

INSTRUCTION MANUAL

7" LCD Digital Desktop
Model# 44353

9" LCD Digital Desktop
Model# 44354

Introduction

Thank you for purchasing a Celestron Labs Digital LCD Desktop Microscope. Your microscope, whether the 7" or 9" model, is a precision optical instrument made from the highest-quality materials. It is designed to give you a lifetime of enjoyment with minimal maintenance.

Before using your microscope, please read these instructions to familiarize yourself with its functions and operations. See the microscope diagram to locate the parts discussed in this manual.

This microscope provides a magnification range of from 120x to 1800x in the 7" screen models and from 152x to 2280x in the 9" screen model. It is well-suited to examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, and more.

Unlike traditional microscopes, the Celestron Labs Digital LCD Desktop Microscope uses an 8MP camera sensor in

lieu of an eyepiece. You will view specimens on a full-color TFT touchscreen measuring 7" or 9". Compared to looking into an eyepiece, using the screen reduces eyestrain and neck strain. The screen also rotates 180° so you can easily share the view with a group.

Onboard firmware allows you to view and capture still images and video of your specimens.

The LCD Desktop Microscopes are packed with traditional high-end microscope features, like fully adjustable mechanical stages, a quad rotating objective turret with four mounted achromatic objectives (4X, 10X, 20X, 40X), and one unmounted (60X) objective. They have coarse and fine focus adjustment knobs, an HDMI cable for outputting full HD video streams to a monitor or projector, a USB-C cable for transferring images to a PC, 10 prepared slides, a four-plug, multi-country AC adapter, a 16 GB MicroSD card, and a dust cover.

Features

- A.** 7" or 9" LCD touchscreen (9" shown)
- B.** Built-in 8MP digital camera
- C.** (1) 4x objective lens, (1) 10x objective lens, (1) 20x objective lens, (1) 40x objective lens
- D.** Full mechanical stage
- E.** Disk diaphragm
- F.** Coarse and fine focus knobs
- G.** Adjustable illuminator

Included Accessories

- H.** AC adapter with four international plugs
- I.** (1) 60x objective lens
- J.** (10) prepared slides
- K.** USB-C cable for transferring images
- L.** HDMI cable for HD-out to monitor/projector
- M.** 16 GB MicroSD card
- N.** Dust cover

Specifications

Type of Microscope:	44353: 7" LCD Digital Desktop	44354: 9" LCD Digital Desktop
Head:	44353: 180° rotatable head - 7" TFT LCD Touchscreen	44354: 180° rotatable head - 9" TFT LCD Touchscreen
Nosepiece:	Quad with click with stop	
Objective(s) Type:	5 achromatic objectives (4 on quad turret, 1 unmounted)	
Objective(s) Standard:	DIN 35 RMS	
Objective(s) Magnification:	4x, 10x, 20x, 40x, 60x (included, not connected)	
Objective(s) Thread:	WJ 4/5" x 1/36"	
Focus Type(s):	Coarse Focus Fine Focus	
Focus Travel (per rotation):	0.2mm to 35.3mm (.007" to 1.39") (est.)	
Magnification Range:	44353: 120x, 300x, 600x, 1200x, 1800x	44354: 152x, 380x, 760x, 1520x, 2280x
Illumination Type:	LEDs	
Illumination Position(s):	Lower, under stage	
Illumination Adjustability Mechanism:	Adjustable wheels at base	
Power Source:	Universal 100 to 240V, 50/60 Hz multi-country plug or 4 AA batteries	
Battery Options:	4 AA, not included	
Stage:	Fully mechanical with metal clips, 88mm x 88mm (3.5" x 3.5") With Filter wheel Red, blue, green 1mm, 3mm, 6mm apertures	
Condenser Type:	Aplanatic	
Condenser N.A.:	N.A. 65	
Stage Diaphragm:	Multiple on wheel	
Eyepiece Type:	CMOS Sensor in lieu of eyepiece	
Eyepiece Magnification:	44353: 30x	44354: 38x
Sensor Size:	8MP 1/2.7"	
Pixel Size:	1.5 µm	
Screen type:	44353: 7" LCD TLT – High-definition touch, 1024 x 600 pixels resolution, 154.21mm (6.07") x 85.92mm (3.38") with 4x digital zoom, rotates 180°	44354: 9" LCD TLT – High-definition touch, 1024 x 600 pixels resolution, 196.60mm (7.74") x 114.15mm (4.49") with 4x digital zoom, rotates 180°
Screen Size:	44353: 7"	44354: 9"
Memory Type:	MicroSD	
Memory Size:	16 GB (card included)	
Port(s) In:	AC Plug (power) Micro SD card slot (supports up to 128 GB) (no card included)	
Port(s) Out:	HDMI USB-C	
AC Adapter:	5-6V DC, 200-2000mA tip positive electrical input	
Still image capture resolution:	(Non-interpolated): 3840 x 2176, 1920 x 1080, 1280 x 736 (Interpolated) 5736 x 3008, 6114 x 3456, 8448 x 4736	
Video resolution:	3840 x 2160, 1920 x 1080, 1280 x 720	
Shutter trigger position:	On LCD Screen	
Software:	Internal Firmware	
Languages Supported	English, French, German, Italian, Spanish, Portuguese, Japanese, Russian, Simplified Chinese, Traditional Chinese	
Microscope Dimensions:	44353: 177.8mm x 177.85mm x 342.9mm (7" x 7" x 13.50")	44354: 388mm x 235mm x 177mm (15.28" x 9.25" x 6.96")
Microscope Weight:	44353: 4.6lbs (2.08kg)	44354: 4.8lbs (2.17 kg)
Included Items:	USB-C cable (data transfer) HDMI cable 10 Prepared Slides 4 plug, multi-country AC adapter Dust cover	

Setup

1. Carefully remove the microscope and other parts from the carton and set them on a table, desk, or other flat surface.
2. Plug the small end from the AC adapter into the socket on the back of the microscope.
3. Choose the correct plug for your country and connect it to the AC adapter. (**Fig. 1**)
4. Plug the AC adapter into your 120-240V power source.

Fig. 1

Using the Microscope

The Celestron Labs Digital LCD Desktop microscope works just like any traditional microscope, except that the image sensor, combined with the LCD screen, replaces the eyepiece.

Power On

Press and hold the orange power button on the back of the LCD screen for 2 to 3 seconds. (**Fig. 2**)

The LCD screen should illuminate, eventually displaying a streaming image.

NOTE: To turn off the microscope, press and hold the power button.

Fig. 2

Using the Illuminator

To turn the illuminator on and adjust its brightness, rotate the illuminator adjustment switch located on the microscope base. (**Fig. 3**)

Fig. 3

Inserting a Specimen Slide

Your microscope ships with a prepared slide kit. Choose one of these slides, or one that you have previously prepared. Your microscope's fully mechanical stage has a slide holder with a spring-activated slide retaining clamp arm. **(Fig. 4)**

1. Open the slide holder mechanism by pivoting the slide retaining clamp arm.
2. Place a glass slide flat on the surface of the stage with the right edge of the slide touching the stationary side of the slide holder.
3. Slowly release the clamp arm to hold the slide in place.

Selecting the Lowest Power Objective

Your microscope is equipped with a quad rotating objective turret. The objectives have powers of 4X, 10X, 20X, and 40X. Rotate the turret so that the 4X objective is in line with the camera sensor and the under-stage illuminator.

Adjusting the Stage

On the left side of the stage, you will see two small black knurled knobs. These knobs are the stage adjustment knobs and control your X and Y movement of the stage. Adjust them while looking at the LCD screen to bring the specimen on the slide into the field of view. **(Fig. 4)**

Focusing

On both sides of the base of your microscope are two different-sized knobs. These are the coarse focus (larger knob) and fine focus (smaller knob). **(Fig. 5)** Start by rotating the coarse focus knob to bring the streaming image of the specimen into a reasonably clear image. Switch to the fine focus and adjust to get your image crisp and clear. Now is a good time to re-adjust the stage position so the specimen on the slide is centered. Now you're ready to begin capturing still images or video—or just continue to observe.

Disk Diaphragm

To refine the beam of light, turn the disk diaphragm and select the one that makes your specimen look best. **(Fig. 6)**

Fig. 4


Fig. 5

Fig. 6

Switching Objectives for More Magnification

To view more details and travel deeper into your specimen, you will have to increase the magnification. Do this by rotating the objective turret until another objective aligns in the light path. **(Fig. 7)** The turret is indexed, so you will feel some tactile feedback when it slots into place. Follow the above focus directions to re-focus on your specimen. If your specimen is filling the complete field of view, you can use your stage adjustment knobs to move along the specimen. The rough total magnification with the five included objectives is based on the magnification of the sensor (44353: 30x; 44354: 38x). Multiply the magnification of the objective by 30 for the total magnification.

Fig.7

44353: 7" LCD Desktop Microscope		
Objective Lens Magnification	Sensor Magnification	Total Magnification
4X	30X	120X
10X	30X	300X
20X	30X	600X
40x	30X	1200X
60X	30X	1800X

44354: 9" LCD Desktop Microscope		
Objective Lens Magnification	Sensor Magnification	Total Magnification
4X	38X	152X
10X	38X	380X
20X	38X	760X
40x	38X	1520X
60X	38X	2280X

Using the Extra 60X Objective

The Celestron Labs Digital LCD Desktop microscope ships with a non-mounted 60X objective. To use it, simply:

1. Unscrew one of the mounted objectives.
2. Align the 60X objective thread with the exposed objective turret thread. **(Fig. 8)**
3. Screw the 60X objective into the turret until snug. Do not over-tighten.

Fig. 8

Using the Touchscreen

The touchscreen should now be active, and the image should be streaming. Familiarize yourself with the main screen interface below. (**Fig. 9**)


Fig. 9

- A. Freeze:** Eliminate camera shake by freezing image before taking a photo.
- B. Settings:** See settings section further on in this manual.
- C. Photo Icon:** Press the photo icon to take a snapshot.
- D. Video Icon:** Press the video icon to begin recording. A timer will appear in the middle/top of the screen, indicating the time elapsed. Press the video icon again to stop recording.
- E. Timelapse Icon:** Alter the timelapse settings in the settings menu and press this icon to begin timelapse.
- F. Gallery Icon:** See gallery section for more information.
- G. Zoom In:** Digitally zoom in to the image. For actual magnification, rotate the object lens
- H. Zoom Out:** Digitally zoom out the image. For actual magnification, rotate the object lens.
- I. Recticle lines Icon:** This places a digital ruler over the image in centimeters. The line color can be changed in the settings menu.

The Settings Menu

1	Exposure	Auto	Line Color	White	9
2	White Balance	Auto	Sound	Off	10
3	Video Resolution	3840x2160	Storage	Micro SD	11
4	Image Resolution	1920x1088	Date Watermark	>	12
5	Black & White	Off	Format	>	13
6	Negative Film	Off	Factory Default	>	14
7	Timelapse	>	Language	English	15
8	Full Screen	Off	Version	>	16

Fig. 10

1. Exposure: Select exposure mode and exposure brightness. (Fig. 11)


Fig. 11

2. White Balance (AWB): Enable one-key white balance and choose automatic or manual mode. (Fig. 12)


Fig. 12

3. Video Resolution: Set the video recording resolution. Click the value to cycle through available resolutions: 3840 x 2160, 1920 x 1080, 1280 x 720

4. Image Resolution: Set the image resolution. Click the value to cycle through available resolutions:

Non-interpolated	Interpolated
3840 x 2176	5736 x 3008
1920 x 1088	6114 x 3456
1280 x 736	8448 x 4736

5. Black & White: Turn the black-and-white photography effect on or off.

6. Negative Film: Turn the negative film photography effect on or off.

7. Timelapse: Choose the duration of your timelapse recording in hours, minutes, and seconds. (Maximum time: 23:59:59.) (Fig. 13)


Fig. 13

8. Full Screen: Enter full screen.

9. Line Color: Change the color of the reticle. (**Fig. 14**)


Fig. 14

10. Sound: Turn the touch sound on or off.

11. Storage: When a MicroSD card is inserted, this screen displays the total memory and available memory. (**Fig. 15**)


Fig. 15

12. Date Watermark:


Fig. 16

13. Format: Format the MicroSD card.

14. Factory Default: Restore the microscope to factory settings.

15. Language: Select language.

16. Version: Display the software version.

The Gallery


Fig. 17

1. Press the gallery icon to view thumbnails of still images and videos you have saved. Click the individual thumbnail to see the full screen.
2. To play a video, tap the screen and then tap the play button.
3. Swipe left/right to scroll through the saved images or use the "Page Back" and "Page Forward" buttons. To get back to the main screen, use the back button.

Using the MicroSD Card

Insert the included MicroSD card or one of your own into the MicroSD card slot if you wish to save images or video (**Fig. 18**)

Fig. 18

Using The HDMI-Out Feature

The Celestron Labs Digital LCD Desktop microscope is equipped with an HDMI cable and an HDMI-OUT port for connecting to a projector or monitor.

1. Locate the HDMI cable and plug one end into the microscope's HDMI-OUT port.
2. Plug the other end into an open HDMI-IN port of the monitor or projector. (**Fig. 19**)
3. Make sure the microscope is powered and turned on.
4. Select the HDMI port as your monitor or projector's source. The view from the microscope will appear on screen.

Fig. 19

Downloading Images To a PC

The Celestron Labs Digital LCD Desktop microscope includes a USB-C cable for transferring images and video to a PC.

1. Make sure the MicroSD card is installed in the microscope.
2. Locate the USB-C cable and plug one end into the microscope's USB-C port. (**Fig. 20**)
3. Plug the other end into an open USB-C port on the PC.
4. The PC will read the microscope (and MicroSD card) as an external drive. Navigate to the drive and open it.
5. Save the images in a location of your choice on your PC's hard drive.

NOTE: You may also use a MicroSD card reader (not included) to transfer images and video.

Fig. 20

Care and Maintenance

Your Celestron Labs Digital Microscope is a precision optical instrument and should be treated with care. Follow these suggestions and your microscope will need very little maintenance throughout its lifetime.

- Store your microscope in a clean, dry place. We suggest you use the included dust cover.
- Never use your microscope where the operation of such devices is restricted. Improper use creates the risk of a serious accident.
- Only use your microscope within the 23° to 120°F (-5° to 50° C) temperature range. Try to limit sudden temperature changes, as this may create moisture in the imager housing.
- Do not attempt to access the internal components of your microscope or modify it in any way. Only authorized technicians should perform maintenance or repairs.
- Keep your microscope away from water and other liquids. Never use it in rain or snow. Moisture creates the danger of fire and electric shock.
- If the exterior of your imager needs cleaning, wipe it with a soft cloth.

Warranty

Your Celestron Labs Digital MicroDirect Microscope has a two-year limited warranty.

Please visit celestron.com for more details.

WARRANTY

www.celestron.com/pages/warranty

FCC NOTICE: This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NEED ASSISTANCE? Contact Celestron Technical Support

celestron.com/pages/technical-support

©2026 Celestron. Celestron and Symbol are trademarks of Celestron, LLC. All rights reserved. • Celestron.com

US: 2835 Columbia Street, Torrance, CA 90503 USA

UK: Unit 2 Transigo, Gables Way, Thatcham RG19 4JZ, United Kingdom

EU Authorized Representative: AR Experts B.V. is Boeingavenue 209, 1119 PD Schiphol-Rijk, The Netherlands, info@ar-experts.eu

This product is designed and intended for use by those 14 years of age and older

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